



Low-height beds

www.muehle-mueller.de

best possible fall prevention • perfect working height

Low-height bed aldena 13.80

A preferably low bed height is especially crucial when considering the increasing number of dementia and Alzheimer's disease and the positioning of patients suffering from these. This applies particularly during sleep as the patients need to be prevented from falling. In order to spare them from the installation of safety bed rails, the company Mühle designed the two nursing beds aldena 13.80 and allegra 20.80. Now the patient is enabled to sleep at lowest height without running the risk of falling out and getting injured, while nursing persons can benefit from the ideal work height of 80 cm!



NEW!

lowest Position 13cm,
highest Position 80cm!

- Enclosed motors, thus concealed mechanics
- Absolute ground clearance after short lift
 - ideal accessibility e.g. for bed lifts
 - easy cleaning conditions below the bed
- Moveable and lockable from any position height
- Comfortable wood casing with integrated safety bars

- Elegantly designed feet and head rests; optional decor
- High safety standard – no crashing or shearing hazards
- Very long life-span (“Made in Germany”)
- Compliance with all standards and directives
- Continuous quality assurance at our own production site



Lowest Position 13cm



Highest Position 80cm



No injury risk when rolling out of bed



excellent fall prevention • perfect work height

Low-height bed allegra 20.80

The nursing bed “allegra” reflects the continuous advancement of our approved nursing bed “rosana”. The low-height bed combines both an extreme low height of 20 cm and a back-friendly work height of 80 cm. Taken as a whole, “allegra” meets all requirements related to modern nursing beds. It distinguishes itself with its elegant appearance, versatility, and the particularly robust design.



NEW!

Height adjustment from 20 up to 80cm!

- Low-height bed with unique adjustability from 20-80cm
- Freely selectable equipment for lying surface
- Modern wood casing with integrated safety side rails
- Optionally available with 2-divided wooden side rails
- Optionally available with bedding rack
- 2-wheels-center lock with 8 double-casters
- Trendelenburg and Anti-Trendelenburg (Feet-down and shock-position)
- Bed extension up to 220 cm
- Very high load bearing capacity (up to 200 kg)
- Stopp- and rollable regardless the height
- High stability and low vibration
- Sustainable economic factor



Lowest Position 20cm



Highest Position 80cm



Max. adjustment of lying surface



Special care beds



Manufacturing

Werner Mühle Metallverarbeitung
Boschstrasse 6
35799 Merenberg

Phone: +49 (0) 6471 - 5065-0
Telefax: +49 (0) 6471 - 5065-13

E-Mail: info@muehle-mueller.de
Web: www.muehle-mueller.de



Sales and distribution

m2 Handels- u. Vertriebs GmbH
Oststrasse 26
51766 Engelskirchen - Runderoth

Phone: +49 (0) 2263 - 9617-0

E-Mail: info@muehle-mueller.de
Web: www.muehle-mueller.de