



mobilia

The Sitting and Rising Bed

mobilia Sitting and Rising Bed

Mobility for patients and easier work for nurses

The mobilia sitting and rising bed in its different versions casa, cura and integra enables the patient to get up from the bed again on his own or with very little support.

Turning lying surface – new mobility

The mobilia nursing bed series is suitable for care at home or in an institution. The effects of many illnesses include physical mobility limitations. Such limitations often affect the patient's ability to move. Turning movements are often impossible, sitting up in bed is difficult or even impossible and will considerably limit quality of life.

The mobilia sitting and rising bed is used for many different disabilities and illnesses, e.g. in case of movement impairment, too-low muscle strength, inability to leave the bed, support at limited stress-capacity of the lumbar spine or joints when sitting due to cardiac insufficiency or for mobilisation, etc.. From the mobilia nursing bed, with all standard functions of a regular nursing bed and an additional sitting position as well as the optional rising aid, patients can better see their environment.

The patient's mobility and transfer are made much easier and perceived as pleasant by the patient. The mobilia sitting and rising beds also give the patient new independence and enable him to rise on his own again. For family performing nursing tasks, the nursing environment and nursing staff, the mobilia series makes things much easier. It prevents occupational diseases in nursing staff, supports family members in nursing and improves the patient's quality of life.

Easy mobility – with the rising bed from the mobilia nursing bed series.

Comfort

- Comfortable lying feeling from
- 2-part mattress
- Rising aid with triangle grip
- Special elastic cover (everything included)
- Arm rests available optionally

Function

- Extreme facilitating of rising and fully motor-powered turning lying surface
- Electrical adjustment of the back, thigh and lower leg rests as well as lifting and rising aid
- Lying surface rotation to either side optional
- Mobilia casa: single wheel braking
- Mobilia cura: single or 2-wheel central braking
- Ground clearance for lifter
- Max. patient weight 100 kg



Design

- Comfortable and beautifully shaped design of the head and foot parts

Areas of use

- Nursing at home, in senior's residences and nursing homes, hospices, rehabilitation facilities and retirement homes, hospices, rehabilitation facilities and retirement homes
- The mobilia **integra** bed insert can be used in personal applications as well as for special uses, e.g. in hotels or cruise ships (for present double or single beds)



mobilia cura

Individuality

- Other wood patterns available on request
- Optionally with steel grid lying surface
- Metal coating in different RAL colours
- Lying surface widths 90 cm or 100 cm
- bed extension to 220 cm possible

Functional Sequence



1. Lying position



2. Head and foot part adjustment



3. Turn of the LS by 45°



4. Turn of the LS by 90°



5. Sitting position



6. Rising position

mobilia cura Shear Lifting Bed

The stable transforming artist

mobilia cura (with its versions mobilia cura, cura e and cura e plus) is one of the two basic versions of our mobilia range. In its basic construction, it is based on a shear-lifting system that particularly represents extreme stability. It also entices with a clear, light-weight design and integrates into any present apartment furnishings.

- Lying surface rotation to either side optional
- Electrically adjustable back, thigh and lower legs rests and height adjustment plus rising aid
- Ground clearance for lifter
- Lying and sitting height smoothly adjustable from 40 - 80 cm (mobilia cura e plus 48 - 80 cm)

Equipment versions



mobilia cura e plus

- Fully electrical lying surface adjustment with rising aid
- Fully electrical turning and rising function



mobilia cura e

- Fully electrical lying surface adjustment
- Fully electrical turning function of the LS



mobilia cura

- Fully electrical lying surface adjustment
- Manual turning function of the LS





mobilia casa Front Lifting Bed

The elegant transforming artist

mobilia casa (with its versions mobilia casa, casa e and casa e plus) is the second basic version of our mobilia range. In this version, the lifting system is integrated into the foot and head parts. This elegant solution offers perfect ground clearance for moving lifter systems under the bed and facilitates cleaning below the bed as well.

- Lying surface rotation to either side optional
- Electrically adjustable back, thigh and lower leg rests, height adjustment and rising aid
- Ground clearance for lifter
- Lying and sitting height smoothly adjustable from 40 - 80 cm (mobilia cura e plus 48 - 80 cm)

Equipment versions



mobilia casa e plus

- Fully electrical lying surface adjustment with rising aid
- Fully electrical turning and rising function



mobilia casa e

- Fully electrical lying surface adjustment
- Fully electrical turning function of the LS



mobilia casa

- Fully electrical lying surface adjustment
- Manual turning function of the LS





mobilia integra Bed Insert

The flexible retrofitting solution for present single or double beds

Physical limitations often happen quite suddenly. A solution must be found quickly to relieve the patient and nurses. For this case, Mühle Müller Pflegebetten offers a mobilia bed version that takes care of such concerns – mobilia integra with all functions of the mobilia system as a pure functional element without bed frame. The mobilia integra can be integrated into any present single or double bed and thus is the perfect solution for anyone who wants to stay in his or her usual home environment.



mobilia integra e plus

- Fully electrical lying surface adjustment with rising aid
- Fully electrical turning and rising function



mobilia integra e

- Fully electrical lying surface adjustment
- Fully electrical turning function of the LS



mobilia integra

- Fully electrical lying surface adjustment
- Manual turning function of the LS





mobilia integra - The perfect integration into
your present bed

mobilia Models

Your individual **mobilia** in a comfortable wood design



napoli



hamilton

All **mobilia** models are optionally available in many high-quality wood and design versions – special wishes will be met as well.

This way, you can design any **mobilia** according to your personal ideas and adjust it to your personal living situation.



hamilton



Production

Werner Mühle Metallverarbeitung
Boschstrasse 6
35799 Merenberg

Phone: +49 (0) 6471 - 5065-0
Telefax: +49 (0) 6471 - 5065-13

E-Mail: info@muehle-mueller.de
Internet: www.muehle-mueller.de



Sales

m2 Handels- u. Vertriebs GmbH
Oststrasse 26
51766 Engelskirchen - Ränderoth

Phone: +49 (0) 2263 - 9617-0

E-Mail: info@muehle-mueller.de
Internet: www.muehle-mueller.de